



ACTGene, Inc.
PO.Box 7011
Piscataway, NJ 08855-7011
Tel: 973-978-6688
Email: sales@actgene.com
Web Site: www.actgene.com

6x HydraGreen™ Safe Loading Buffer

Product contents:

P/N: ACT-IDMG05

1x1ml or 1x20ul (Sample), **1ml 6x Loading Buffer sufficient for 500-1000 sample loadings**

Introduction:

6x HydraGreen™ Safe Loading Buffer is a non-carcinogenic safe gel loading buffer containing density agents, tracking dyes, and HydraGreen™ Safe DNA Dye. The 6X loading buffer is added to samples in place of gel loading buffer, and eliminates the need to add DNA dye to the agarose gel during casting or staining buffer after gel running.

HydraGreen™ Safe DNA Dye (included in 6x HydraGreen™ Safe Loading Buffer) has fluorescence excitation maxima at 295 nm and 490 nm. The fluorescence emission maxima is similar to EtBr when bound to DNA – at 530 nm.

Protocol:

Note: It is highly recommended agarose gel to be ≤ 0.5 cm in thickness. TAE buffer system is preferred for higher sensitivity.

Prepare agarose gel according to standard protocol. Do NOT add any dye to the agarose or buffer.

Mix 6x HydraGreen™ Safe Loading Buffer in a **1:5 ratio** with your DNA samples. For example, for 10 μ l of DNA sample, add 2 μ l of 6x HydraGreen™ Safe Loading Buffer or for 5 μ l of DNA sample, add 1 μ l of 6x HydraGreen™ Safe Loading Buffer.

Mix thoroughly before loading.

DETECTION :

Detect bands under UV illumination (yellow or green gelatin- or cellophane filters is recommended for clearer bands) or non-UV LED illuminators such as Blue Light LED illumination.

Shipping and Storage Conditions :

Shipped in RT

Store in 4 °C (protect from light)

For Research Use.