


ACTGene™ DNA Marker

Innovative Nucleic Acid Electrophoresis Workflow Essentials

DNA Markers


- 100 bp, 1 kb, & 1 kb Plus marker choices
- Fragments with bolded intensity for easy viewing
- Ready-to-use format
- Idea for agarose gel electrophoresis, Northern, & Southern applications



1.5% Agarose Gel 1xTAE
(100-1,500 bp) 11 Fragments



1% Agarose Gel 1xTAE
(250-10,000 bp) 13 Fragments



1.5% Agarose Gel 1xTAE
(100-10,000 bp) 19 Fragments

Order Information

ACT-IDMWD100 ACTGene™ DNA Marker, 100 bp Ladder, 11 Fragments, 500 ul (sufficient for ~100 loads)

ACT-IDMWD1 ACTGene™ DNA Marker, 1 kb Ladder, 13 Fragments, 500 ul (sufficient for ~100 loads)

ACT-IDMWD2 ACTGene™ DNA Marker, 1 kb Plus Ladder, 19 Fragments, 500 ul (sufficient for ~100 loads)

ACTGene, Inc.

Tel: (973) 978-7888 | Fax: (732)469-0333

E-mail:sales@actgene.com