



ACTGene, Inc.
P.O. Box 7011
Piscataway, NJ 08855
Tel: 973-978-7888
Fax: 732-469-0333
Email: sales@actgene.com
Web Site: www.actgene.com

Deoxynucleoside Triphosphate (dNTP)

Catalog Numbers:	<u>Cat Number</u>	<u>Deoxynucleotide</u>
	2013MIX	dNTP Mix
	2013SET	dNTP Set
	2013A	dATP
	2013G	dGTP
	2013C	dCTP
	2013T	dTTP

Product Description: dNTP Mix (p/n: 2013MIX) is premix solution containing 25mM of each, supplied in one tube, with 100mM total conc. at 1ml.
dNTP Set (p/n: 2013SET) is a set of dNTP consisting of 4 separate tubes (4x1ml) each with 100mM concentration of dATP, dGTP, dCTP, and dTTP.
dATP (p/n: 2013A), dGTP (p/n: 2013G), dCTP (p/n: 2013C), dTTP (p/n: 2013T) is supplied in individual tube containing each at 100mM conc. at 1ml.
(Customizable concentration and packaging available upon request).

Storage Condition: -20°C or -70°C constant temperature freezer.

Purity: dNTP and dNTP mix are >99% chromatographically pure. Free of DNase, RNase, Protease, or Phosphatase activities.

Application: For direct in-vitro synthesis of DNA