



Product Specification

Speed	10-40 rpm
Environmental temperature	4-60 °C
Max load	0.8kg (include rotisseries)
Time for acceleration /deceleration	16 sec to and from highest speed
Centrifugal medium density	1.2g/ml
Weight	1.2 kg
Dimension	22.5x15.4x19.5 cm(w x d x h)
2 Modes Selection	Variable speed & Oscillating

Product Highlights

- ➔ Digitally controlled variable speed and oscillating modes
- ➔ Oscillating mode features 40° rotation, pause, then back and forth five times before another 40° to repeat the process.
- ➔ Accommodating 0.5-50ml tubes
- ➔ Digital key pad for easy speed adjustment 10-40 rpm
- ➔ Small footprint, modern compact design
- ➔ Super quiet
- ➔ Can use 4 rotisseries simultaneously
- ➔ Can be used in incubator 60° or cold room 4°
- ➔ Optional Magnetic platform for Western Blot Application

Ordering Information

- ACT-1040** Piggy™ Tube Rotator, including rotisserie for 42 x 0.5mL~0.8mL tubes, 2/pk, 50 x 1.5mL~2.0mL tubes, 2/pk, rotisserie for 14 x 5mL~7mL tubes, 2/pk, rotisserie for 14 x 10mL~15mL tubes, 2/pk, rotisserie for 6 x 50mL tubs, 1/pk
- ACT-1040-6** Magnetic platform, accessory for Piggy™ Tube Rotator, includes 4 metal mounting strips

ACTGene, Inc.

P.O. Box 7011

Piscataway, NJ 08855

Tel: (973) 978-7888 | Fax: (732)469-0333

E-mail: sales@actgene.com