

## ACTGene™ Protein Marker 10-245 kD, 12 Protein 3 Color

### Product contents:

P/N: ACT-IDMW24

1x500µl

Sufficient for approximately 150 loads

---

### Introduction:

ACTGene™ Protein Marker is twelve protein three color pre-stained marker covering wide range of 10-245kD protein sizes. Conveniently delivered in ready-to-use format and ideal for SDS-PAGE, Western Blotting, and related protein applications.

---

### Protocol:

Swirl tube well.


Add ~4µl per load onto protein gel prior to electrophoresis.

---

### Shipping and Storage Conditions:

Shipped in RT

Store in in 4°C for short term and -20 °C for long term



Tris-Glycine 4-20%  
(10-245kD) 12 Protein 3 Color