



# AgilePipettor

## Single Channel Adjustable Volume Micropipettors



ACTGene goal is to offer high-tech products, at the best quality, giving added value and innovative features to the user, all at an equitable price. AgilePipettor is the volumetric instrument designed to measure and transfer liquids precisely and safely. It can measure and transfer volumes from 1 $\mu$ l to 1000 $\mu$ l.

### Features :

- Fatigue free working with convenient finger rest, AgilePipettor features significantly lower pipetting forces than the average pipetting forces in the market.
- Made with high-density, chemical-resistant polymer for better protection against chemical and physical corrosion.
- Volume display is always visible while pipetting.
- Locking mechanism preventing accidental volume changes during pipetting.
- Light & smooth plunger action ensures precise aspiration.

### Specifications :

Cat. No.	Volume Range(ul)	Increment (ul)
AMP-10	1 ~ 10	0.02
AMP-20	2 ~ 20	0.02
AMP-200	20 ~ 200	0.2
AMP-1000	100 ~ 1000	2.0

