



## Product Specification

Speed	12,000 rpm max.
RCF	7,500 x g
Rotor (standard)	1) Round rotor: 6-place for 1.5mL / 2.0mL tubes 2) 6-tube adapter for 0.2mL or 6-tube adapter for 0.5mL tubes 3) Strip rotor: 16-place for 16 x 0.2mL tubes or 2 x 0.2mL 8-tube PCR strips
Rotor (optional)	Optional Round rotor: 6 x 1.5mL/2.0mL + 6 x 0.5mL + 6 x 0.2mL
Weight	3.5 kg
Electrical	85~245V AC, 50Hz,
Environment	Can be operated in 50°C and 4°C environment

## Product Highlights

- ➔ Maximum speed: 12,000rpm(7,500xg)
- ➔ Small footprint, modern compact design
- ➔ Supplied complete with rotor and adaptors
- ➔ Accelerates and stops quickly
- ➔ Super quiet
- ➔ Build in safety - Rotors stop once lid release button is pushed; rubber feet for vibration-free operation
- ➔ Ideal for quick spins

## Ordering Information

ACT-12000 ACTGene™ Microcentrifuge, 12,000RPM



Standard Rotors and Adaptors



Optional Rotor

**ACTGene, Inc.**

P.O. Box 7011

Piscataway, NJ 08855

Tel: (973) 978-7888 | Fax: (732)469-0333

E-mail: sales@actgene.com