

ACTGene™ DNA Marker 1 kb Ladder, 13 Fragments

Product contents:

P/N: ACT-IDMWD1

1x500µl

Sufficient for approximately 100 loads

Introduction:

ACTGene™ DNA Marker is thirteen DNA fragment ladder covering wide range of 250-10,000 bp sizes. Conveniently delivered in ready-to-use format and ideal for agarose gel electrophoresis, Northern, Southern, and related DNA analysis applications.

Protocol:


Swirl tube well.

Add ~5µl per load onto agarose gel prior to electrophoresis.

Shipping and Storage Conditions:

Shipped in RT

Store in in 4°C for short term and -20 °C for long term



1% Agarose Gel 1xTAE
(250-10,000 bp) 13 Fragments