

AgilePipettor™ Electronic Single Channel Pipettor



Product Highlights

- ➔ Forward Pipetting
- ➔ Sequential Dispensing
- ➔ Reverse Pipetting
- ➔ Multi Dispensing
- ➔ Manual Pipetting
- ➔ Multi Aspiration
- ➔ Mixture
- ➔ Pipetting + Mixture
- ➔ Autoclavable shaft
- ➔ **Clear LCD screen display**
- ➔ **Ergonomically designed**
- ➔ **Compatible with universal fitting tips**

Product Highlights

Part Number	Volume range (µl)	Test volume (µl)	Accuracy (±)	Precision (±)	Tip type (µl)
ACT-ESCP-10 (.5-10ul)	0.5-10	0.5	0.020	0.010	10
		1	0.025	0.020	
		5	0.075	0.040	
		10	0.100	0.040	
ACT-ESCP-100 (5-100ul)	5-100	5	0.10	0.10	200
		10	0.20	0.10	
		50	0.50	0.15	
		100	0.80	0.20	
ACT-ESCP-300 (20-300ul)	20-300	20	0.50	0.14	200
		30	0.75	0.21	
		150	1.50	0.45	300
		300	1.80	0.60	
ACT-ESCP-1200 (50-1200ul)	50-1200	100	3.00	0.60	1000
		120	3.60	0.72	
		600	6.00	1.2	2000
		1200	7.20	2.4	

ACTGene, Inc.

P.O. Box 7011

Piscataway, NJ 08855

Tel: (973) 978-6688

E-mail: sales@actgene.com

www.actgene.com