



Palm™ Dry Bath, Portable with optional car charger



Product Specification

Temperature / Accuracy	-10 to 100 °C / $\leq \pm 0.5^\circ\text{C}$ at 100°C
Time range	0 – 999 min.
Cooling time	≤ 15 min. (from RT to RT-15 °C if RT < 30°C) ≤ 30 min. (from RT to RT- 20 °C if RT < 25°C)
Thermo Lid	YES
Weight	1.4 kg
Dimension	11x15.6x100 cm
Max heating power	60W
Power supply	AC 110-220, optional car charger

Product Highlights

- Palm size 11x15.6x10cm
- Light weight of only 1.4kg
- **Fully portable with optional car charger**
- Heating and Cooling functions in one system
- For sample size up to 2ml
- Lid heating function to prevent sample evaporation
- Over-heating protection safety mechanism
- Build-in temperature calibration function
- Various blocks available

Ordering Information

ACT-PALM-24	Palm™ Dry Bath, Portable, -10 to 100 °C 1.5kg (main unit only, blocks optional)
ACT-PALM-24A	40 x 0.2mL block / 8 x 0.2mL strip
ACT-PALM-24B	24 x 0.5mL block
ACT-PALM-24C	15 x 1.5/2.0mL block
ACT-PALM-24D	8-12.5 x 12.5mm cuvette block
ACT-PALM-24E	4 x 15mL block (lid can't be closed)
ACT-PALM-24F	2 x 50mL block (lid can't be closed)
ACT-PALM-24G	Car charger

ACTGene, Inc.
P.O. Box 7011
Piscataway, NJ 08855
Tel: (973) 978-6688
E-mail: sales@actgene.com